

GIROPHARMA TITAN CARROUSEL STORER GG



OPTIMIZE

your storage space and increase your productivity




SECURE

your product removal and access



SIMPLIFY

your management and avoidout of stock



Up to 90% less space occupied



Accident prevention



Controlled environment (humidity, temperature, ISO norms and fire resistance)



Reduction of expiries



Badge and password identification



Over 99% availability rate



Less picking errors



Standard or custom project-based solutions



Traceability of the products



Management interface possible with your ERP



Management clean room environnement - CH de MONTLÉGIA

Implementation of carousel storers with ISO standards.

Objective: Management of the stock of drugs and SDM within an operating theater.



Control of the drugs temperature - CH LA PITIÉ SALPÉTRIÈRE

Implementation of a controlled temperature between 2 and 6°.

Objective: Management of stocks of drugs dedicated to chemotherapy.




Total safety of blood bag movements - EFS


Secure access by badges and computerized management of the blood bag stock.

Objective: Identification and total security of blood bags within of the French Blood Donors Organisation.

Dynamism, expertise, analysis, innovation



Key figures



Contact us!



ELECTROCLASS
12 avenue Gutenberg,
PA Gustave Eiffel
77600 Bussy-Saint-Georges

t. +33 01 80 66 99 90
f. +33 01 64 76 63 72
info@electroclass.com
www.electroclass.com

Simplified joint-stock company (French SAS) with a registered capital of 546 100 €
Meaux Trade and Companies Register (French RCS) 2008B02121

APE business identifier: 4666Z
SIRET number: 316 654 441 00096

Intracommunity VAT number:
FR 43 316 654 441




TECHNICAL SPECIFICATIONS

 Load-bearing chain taking up to 500 kgs effective load per shelf


 Frequency-drive




Clear, readable 5.7' color touch screen



EU N°15095+A1 standard, class 4 safety light curtains



Area indicators with depth indicator



AC motor

 Work station LED lighting

GIROPHARMA RANGE

	G 40 / G 50	HD40	G 60	G 80	G 160
Height (m)	From 2.4 to 6	From 2.4 to 6	From 2.5 to 7	From 2.5 to 7.5	From 2.7 to 10
Overwall depth (mm)	1.180	1.560	1.725	1.596	1.830
Lookup table depth (mm)	300	300	300	300	400
Effective shelf depth (mm)	From 340 to 400	From 520 to 620	From 520 to 620	From 420 to 520	From 520 to 620
Overwall width (mm)	Useful + 400	Useful + 400	Useful + 610	Useful + 758	Useful + 860
Effective width (mm)	From 1820 to 2880	From 1820 to 2880	From 1950 to 2950	From 2480 to 2900	From 1400 to 3500
Max. load per shelf (kg)	150	150	250	300	500

G-STOCK PHARMA WMS software

G-STOCK PHARMA, published by TEMREX, is best-in-class WMS for managing your storers in a demand-responsive way as well as:

- > managing stock and spaces
- > managing items, packaging types, etc.
- > tracking all movements
- > managing batches, expiry dates, serial numbers, secure access, WiFi, etc.


A 100% TEMREX service team

- > A national service network and a wide availability of parts to ensure a quick response time
- > A hotline for remote diagnosis and problem solving in man cases
- > Maintenance contracts adapted to ensure maximum availability

100% TEMREX solutions

Drawing on its experience and rigor, our team of designers today harmoniously combines technique and user-friendliness to offer you carousel storers dedicated to enhancing your performance.

Through this approach our high-availability, connected solutions improve the productivity of tomorrow's industry.



Ventilation, temperature, fire protection

GIROPHARMA storers are frequently equipped with ventilation systems (CMV), with systems for ensuring temperature and humidity remain constant as well as a cooling unit, etc.



Managing products of heterogeneous volumes

To make the handling of medical devices or boxes of medicines easier, GIROPHARMA storer are easily modular. Arrangements (molded containers, dividable bins) will allow you to configure the storage unit according to your needs.

